

Molina Healthcare of South Carolina Provider Update



Important Information — Please Read

PREFERRED DRUG LIST UPDATE / PRODUCT DISCONTINUATION ALERT

Dear Provider:

On April 11, 2016, the manufacturer Trividia Health informed Molina Healthcare that it is discontinuing its diabetic testing products TRUEtest[®] test strips and TRUEresult[®] and TRUE2go[®] blood glucose systems in the United States. The TRUEtest[®] test strips are GDH-PQQ enzyme-based test strips. For a very small percentage of diabetics, the GDH-PQQ enzyme can interfere with certain drug products or therapies that contain non-glucose sugars.

TRUEtest[®] test strips and TRUEresult[®] blood glucose system have been the exclusive, preferred brands of diabetic testing supplies for Molina Healthcare of South Carolina. We along with Trividia Health are committed to providing the best possible patient and provider support during this time to ensure a smooth conversion from these discontinued products.

Here are some important points to note:

- Trividia is supporting a free upgrade for patients to one of their newer brands, TRUE METRIX[®] or TRUE METRIX[®] AIR blood glucose system.
 - TRUE METRIX[®] AIR has the latest technology. This meter is now Molina's preferred product for this conversion for our members, effective immediately.
 - Enclosed is a product information guide that compares the features of the new TRUE METRIX[®] AIR blood glucose system to those of the discontinued TRUEresult[®] meter.
- Molina Healthcare is working closely with its pharmacy benefit manager vendor, CVS Health, and Trividia Health to assist with appropriate and timely meter and letter mailings for impacted members. *No action is required on your part to convert your patients to the TRUE METRIX[®] AIR meter.*
- This is NOT a product recall. This is a manufacturer's voluntary product discontinuation.
- To assist you when prescribing testing supplies for your diabetic patients whom are Molina Healthcare members, below is a table of the discontinued product NDCs and the NDCs for the products to which Molina is converting.

PRODUCT LIST OF DISCONTINUED NDCs		
PRODUCT	NDC	NOTES
TRUEresult [®] Blood Glucose System	56151-1240-01	Placed an NDC block on this product
TRUE2go [®] Blood Glucose System	56151-1340-01	
TRUEtest [®] Test Strips (25 count)	56151-1030-25	
TRUEtest [®] Test Strips (50 count)	56151-1030-50	Will process as preferred while supplies last
TRUEtest [®] Test Strips (100 count)	56151-1030-01	Will process as preferred while supplies last
NEW PREFERRED NDCs for MOLINA HEALTHCARE		
PRODUCT	NDC	NOTES
TRUE METRIX [®] Blood Glucose System	56151-1470-02	
TRUE METRIX [®] AIR Blood Glucose System	56151-1490-02	Added to the Molina PDL as preferred effective immediately
TRUE METRIX [®] Test Strips (50 count)	56151-1460-04	
TRUE METRIX [®] Test Strips (100 count)	56151-1460-01	




Molina Healthcare of South Carolina Provider Update

Important Information — Please Read



Here is a timeline of events as it pertains to this product discontinuation:



IMPORTANT PHONE NUMBERS	
Molina Healthcare of South Carolina Pharmacy Prior Authorization:	TEL: 1-855-237-6178
Hours of operation: 8:00 a.m. to 6:00 p.m. EST Monday – Friday	FAX: 1-855-571-3011

If you have questions regarding these changes, please call Molina Healthcare of South Carolina at 1-855-237-6178.

Thank you for your partnership in providing quality care for Molina Healthcare of South Carolina members.

Sincerely,
Pharmacy Services Department
Molina Healthcare of South Carolina

Enclosure





Product Features	TRUE METRIX® AIR Self-Monitoring Blood Glucose System	TRUEresult® Blood Glucose Meter
Coding Required	No	No
Sample Size	0.5µL	0.5µL
Test Speed	As fast as 4 sec	As fast as 4 sec
Ketone Test Reminder	Yes	Yes
Weekday Display	No	No
Averaging	7/14/30/60/90	7/14/30
Event Tags	Yes (6 sets)	No
Downloadable Results	Yes	Yes
Strip Release	Yes	Yes
Testing Reminder Alarms	Yes (4 sets)	Yes (4 sets)
BG Range	20-600	20-600
Memory	1000	500
Strips Used	TRUE METRIX® Self Monitoring Glucose Test Strips	TRUEtest Strips® Blood Glucose Test Strips
Hematocrit Range	20%-70%	20%-60%
Temp Range (system ops)	50°-104°F	50°-104°F
Humidity	10-90%	10-90%
Chemistry	GDH-FAD	GDH-PQQ
Blood Sampling	Capillary whole blood	Capillary whole blood
Neonate Claim	No	No
Venous Claim	Yes	Yes
Arterial Claim	No	No
Wireless Technology	Bluetooth Smart	NA
App	TRUE MANAGER AIR App	NA
App Operating System Availability	Android and iOS	NA

TRUE METRIX AIR



FEATURES:

- No Coding
- Tiny 0.5 microliter sample size
- As fast as 4 seconds
- Stores 1,000 results with time/date
- Data management
- Download capabilities

with



Featuring



Includes Event Tagging



BEFORE MEAL
Test was taken just before a meal



AFTER MEAL
Test was taken 2 hours after the start of a meal



EXERCISE
Test was taken during or after exercise



SICK
Test was taken when sick



MEDICATION
Medication taken may have affected test result



OTHER
Any other reason that the test is unique